

Magnetic stirrer MS-01

MS-01 is a powerful, efficient, stainless steel magnetic stirrer.



FEATURES

Strong: This motor-driven unit provides an extra strong field for up to 4 flat-bottom lab flasks at one time.

Versatile: Specially designed to be compatible with the TW-2.03 water thermostat, the combo offers highly accurate temperature control and effective mixing simultaneously.

Durable: The magnetic stirrer's stainless steel enclosure is highly resistible, as well as easily cleaned and sterilized.

Convenient: Its simple and noiseless operation make the MS-01 an obvious choice for any laboratory.

MAGNETIC STIRRER MS01 SPECIFICATIONS

Speed range	RPM	1 - 2200
Liquid volume range for one mixing place	L	from 0.2 up to 2L
Power adapter output	$\ominus \text{---} \oplus$	DC 12V, 1.5A
Dimensions (length x width x height)	mm	260 x 260 x 40
Weight	kg	2.1